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24.09.03 CB/GP

Operating Manual Mode d'emploi

Peristaltikpumpe

BV - GES

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Peristaltic pump

BV - GES

multi channel, with interchangeable tube-bed sets
ISM 899

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Français

The ISMATEC logo features a stylized graphic of several overlapping, rounded rectangular shapes in shades of gray above the word 'ISMATEC' in a bold, sans-serif font, followed by a registered trademark symbol (®).

ISMATEC®

The IDEX Corporation logo features the word 'IDEX' in a bold, italicized, sans-serif font above the words 'IDEX CORPORATION' in a smaller, standard sans-serif font.

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Contents

Safety precautions	3–5
Warranty terms	6
Product	7
Operating panel	8
Rear panel	8
Mains voltage	9
Voltage setting	10
Changing the fuses	10
Tube-bed sets	11–13
Inserting the tubing	14
Starting the pump	15–16
Pumping against back-pressure	17
Overload protection	18
When the pump is not in use	19
Running-in period for tubing	20
Tubing life	20
Analog interface	21–22
Settings of switz S1	23
Accessories	24
Flow-rates	25
Maintenance	26
Repairs	26
Spare-parts	26
Technical specifications	27

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pump-head). Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used. The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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✉ Warranty

For all parts manufactured by ISMATEC®

from date of delivery: 2 years

All other parts,

excluding consumables: 1 year

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC SA the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in case of

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing

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Product

Package contents

- Peristaltic pump BV-GES (drive)
- Tube-bed sets SB-2V and/or SB-3V as ordered
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Operating panel

- 1 Mains switch
- 2 Speed selector
999 = max. speed, adjustable
in 0.1% steps
- 3 Rotation direction
- 4 MIN / MAX key for speed
(ideal for filling or emptying
the tubes)
- 5 Run / Stop

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
⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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
Voltage setting 115V / 230V and changing the fuses

- Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by opening
it with a small screw-driver (size
0).
 - 2 Take out the voltage selector
plate. Turn it and re-insert it into
the fuse-holder so that the
required voltage rating is facing
the window of the fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x 1.25 A (slow-blow)
115 V_{AC}: 2 x 2.50 A (slow-blow)
⚠ Use always 2 slow-blow fuses
of the same type complying with
the local mains voltage.
 - 4 Shut the fuse-holder. The voltage
rating is visible in the window.

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

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Tube-bed sets

SB 2V ISM 010 (2 channels)
 Type of tube: Standard
 Tube i.d.: 3.2 to 8.0 mm
 Flow-rates: 2.6 – 1110 ml/min

SB 3V ISM 011 (3 channels)
 Type of tube: Standard
 Tube i.d.: 0.8 to 4.8 mm
 Flow-rates: 0.22 – 530 ml/min

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Examples of different tube-bed combinations (with 4 to 6 channels)

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
Mounting the tube-bed sets SB 2V and SB 3V

On each side of the pump-head one of the 2 tube-bed sets SB 2V or SB 3V can be mounted.

1 Insert the 2 or 3 tube-beds.

- a) Thick bolt for tube-beds
- b) Thin bolt for pressure levers

2 Secure the tube-beds with the large flat fillister head screw.

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Mounting the tube-bed sets SB 2V and SB 3V

3 Mount the 2 or 3 pressure levers

➔ Do not forget the spacer rings!

➔ Each spring-loaded bolt of the 2 or 3 pressure levers must be pointing towards the vertical center of the corresponding tube-bed!

4 Secure the pressure levers with the small flat fillister head screw

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
Inserting the tubing

➔ Switch the pump off
Open the pump-head protection shield. Open the tube-bed – turn the pressure lever in the clockwise direction

1 Push both ends of the tubing through the tubing retainers

2 Stretch the tubing alongside the roller-head and fix it with the knurled screw

➔ Tighten the screw only as much as is necessary in order to prevent the tubing from moving at the required rotation speed

➔ Press the tube-bed towards the tubing and close the pressure lever by turning it to the left

➔ Close the pump-head protection

3 Adjust the pressure onto the tubing by turning the screw whilst the pump is running – for each pump channel individually – see Fig. 3

4 Tube-beds with completely inserted tubing

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
Starting the pump

- Check the voltage rating in the window (A) of the fuse-holder (on rear panel).
If necessary, change the rating as indicated on page 10.
- Mount the tube-beds according to the instructions given on pages 12 to 13.
- Insert the tubing according to the instructions given on pages 14.
- Connect the tubes to the system.
- Connect the pump to the mains.

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Starting the pump

- 1 Press in the mains switch.
- 2 Set the required speed on the 3-digit speed selector.
(999 = 240 rpm)
➔ The rotation speed can be adjusted while the pump is running.
- 3 Choose the rotation direction
➔ Never change the rotation direction while the pump is running.
- 4 The MAX/SLOW- button allows rapid filling (or emptying) of the tubes while the pump is running.
- 5 Start the pump by switching the RUN/STOP button.

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Pumping against pressure

The BV-GES can be used for continuous duty at a max. differential pressure of about 1.5 bar (21 psi). To a great extent this value depends on the max. operation pressure of the tubing used. Depending on the tubing material and the relation between the bore size and the wall thickness, the max. differential pressure may be lower or higher.

Apart from the tubing material the max. differential pressure is also determined and limited by the performance (torque) of the pump motor as well as the design of the pump-head. In addition, short-time operations may generally allow higher differential pressures than long-term duties.

In case of any doubts please contact your ISMATEC® agent.

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Overcurrent protector

The tubing pump BV-GES features an overload protector.

When the overcurrent protector is activated, the pump is stopped. In case of an overload, pin 12 of the analog interface (open-collector) is connected to the ground (pin 1) (see page 22).

Whenever this situation occurs, **the pump must be switched off immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again)

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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When the pump is not in use

When the pump is idle, we recommend you to release the tubing from pressure.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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⚠ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used. Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Analog interface

Pin 1, (GND) Ground

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND)

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

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
Analoginterface

Pin 7, +36 V_{DC}

About +36 V_{DC} (unstabilised) are available the pump motor (max. load 1 A).

Pin 9, speed OUT

The default setting is 0–10 V_{DC}, proportionally to the motor speed 2.4–240 rpm. Alternatively a frequency range from 0 to 12 kHz is available. The input range can be selected via a dip-switch inside the pump (see page 23).

Pin 10, +5V_{DC}

This pin is connected to the stabilised +5V_{DC} voltage supply of the pump.

Pin 12, overload

In case of an overload, this pin (open-collector) is connected to the ground.

➔ Switch off the mains voltage and wait 2 minutes for cooling down.

Pin 13, internal speed

If pin 13 is connected with pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at pin 5.

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
Broche 13, vitesse internal (internal speed)

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Pin 9 Ausgang/Output/Sortie
 A = 0–10 V_{DC} (Standard)
 B = 0–12 kHz

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Pins		DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN Entrée vitesse	0–5V	100kΩ	ON*	OFF*	OFF*	OFF*	ON*
	0–10V	20kΩ	OFF	ON	OFF	OFF	ON
	0–20mA	250Ω	OFF	OFF	ON	OFF	ON
	4–20mA	250Ω	OFF	OFF	OFF	ON	OFF

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➔ **Make sure that the pump is disconnected from the mains supply.**

⚠ The instrument should only be opened by a qualified technician!
 Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Foot switch

Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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⚠ Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®

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The indications in the table are only approximate values and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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Maintenance

Provided the BV-GES drive is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective BV-GES drive to your ISMATEC® representative.

Please give information on defect and date of purchase.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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△ Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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IP 30

Operating conditions

Temperature +41 to +104°F
Rel. humidity max. 80%
– not condensing, at normal laboratory conditions

Dimensions/Weight (Drive)

DxWxH 7.5x6.9x9.85 inches
Weight 14.3 lb

CE-compatibility proved according to: EN 61326-1/EN 61010-1

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
Ask for the latest
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH