ISMATEC®

IDEX
HEALTH & SCIENCE

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15.03.07/CB/GP

Operating Manual Mode d'emploi

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Ecofluidic

1 Kanal, 3 Rollen
ISM 1200

Deutsch

Tubing pump

Ecofluidic

1 channel, 3 rollers
ISM 1200

English

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Ecofluidic

1 canal, 3 galets
ISM 1200

Français

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Contents

Safety precautions	4
Warranty terms	7
Product	8
Rear panel	9
Mains voltage	10
Voltage setting	11
Changing the fuses	11
Starting the pump	12
Control elements	12
Insert the tubing	13
Operation	14
Running-in period for tubing	15
Tubing life	15
Pumping against pressure	16
When the pump is not in use	16

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Contents

Overload handling	17
Status Indication	19
Flow rates	22
Maintenance	23
Repairs	23
Spare-parts	23
Disposal	23
Cleaning	24
Technical specifications	24

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.

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Safety precautions

- The pumphead must be opened only while the pump is at rest. If it is opened on a running pump, the safety stop is activated. This protection relies on a contactless sensor that is located above the closure lever. **The safety stop must not be deactivated by manipulations on the pump!**
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.

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Safety precautions

- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping aggressive media, we recommend to place the pump elevated in a vat with outlet.
- Only new fuses must be used. They have to comply with the specifications given on page 10. In case of failure, always both fuses must be exchanged.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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✉ Warranty

**For all parts manufactured by ISMATEC®
from date of
delivery: 2 years**

**In case of any queries, please contact
your ISMATEC® representative.**

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packing or in a packet of equal quality. The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Package contents

- Tubing pump Ecofluidic
1 channel, 3 pump rollers
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 Mains socket
- 2 Mains switch
on/off (start/stop)
- 3 Opening for screw-driver
- 4 Fuse-holder containing 2 fuses
(change always both fuses)
- 5 Voltage selector

⚠ Ventilation slots (a)

For sufficient ventilation keep the slots at the site and the base slots clear.

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⚠ Socket/Power cord

Use exclusively the originally supplied power cord.

The socket must be earthed (protective conductor contact).

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

- Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by opening it
with a small screw-driver (size 0).
 - 2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse-
holder so that the required voltage
rating is facing the window of the
fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x 1.6 A (slow-blow)
115 V_{AC}: 2 x 1.6 A (slow-blow)
⚠ Use always 2 slow-blow fuses of
the same type complying with the
local mains voltage.
 - 4 Shut the fuse-holder. The voltage
rating is visible in the window.

⚠ Before starting-up

Check if the voltage setting visible in the
window of the fuse-holder complies with your
local mains voltage.

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
Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel).
If necessary, change the rating as indicated on page 11.
- Insert the pump tubing (see page 13).
- Connect pump to mains and switch it on at the back panel.

Control elements

1. Control knob
2. Status LED
3. Stand-by position (LED = yellow)
4. Speed in ccw direction
5. Speed in cw direction
6. Start / stop by pressing the control knob (LED = green)

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
Insert the tubing

→ Stop the pump. The LED color is yellow (stand-by)

⚠ If the LED flashes yellow, the lid had been opened while the pump was running. Consequently, the safety monitor has stopped the pump. Press the control knob once to reset it to 'stand-by'.

- 1 Push the lid down at the upper left indicator point and open the lid lock (Fig. 1). Swing the lid open. In the beginning, you will feel resistance from the snap-in point.
- 2 Untighten the knurled head screw completely. Lift the tubing clamp and turn it upright (Fig. 2). Rotate the rollerhead by hand in 12 o'clock position.
- 3 Insert tubing without twisting it and let it hang down. Close the tubing clamp but do not tighten it completely yet. Gently pull the tubing and fix the tubing clamp now (Fig. 3).
- 4 Close the lid beyond its snap-in point (Fig. 4) and secure it with the lid lock. If the lid gets blocked at half-way by the knurled head screw, the tubing clamp has not been tightened properly.

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Operation

Press the operation knob

Start, LED changes from yellow to green

Turn operation knob clockwise

Pumphead rotates cw with increasing speed

Press operation knob again

Stop, LED changes from green to yellow as soon as the pump has stopped rotating (ca. 1 - 2 sec)

Press operation knob once again

Start, the pump runs at the same speed as before

Turn operation knob while pump is idle

Speed preset for next start

⚠ Under high load conditions the LED will possibly flash yellow-green right after the start for a short time. When reaching its set speed, it changes to steady green.

⚠ When using new tubing it can happen that the pump initially doesn't move the medium, depending on the tubing (hardness and diameter). If this is the case, we recommend to make the inside of the tubing wet and run the pump with inserted tubing for 15 to 30 min.

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⚠ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used: Chemicals, rotation speed, differential pressure, temperature, viscosity, etc.

General information on the service-life is stated in our detailed tubing documentation (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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Pumping against pressure

The Ecofluidic can be used for continuous duty at a differential pressure of max. 1.5 bar (smaller tube sizes allow higher pressures)

In case of any doubts please contact your ISMATEC® agent.

If the pump drive is subjected to excessive load, overload protection comes into action, see page 17.

When the pump is not in use

When the pump is idle, we recommend to release the tubing from being squeezed. This helps to protect the tubing from unnecessary strain, maintains its suction capability and prolongs its service-life.

Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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Overload handling

The pump Ecofluidic has an overload protection with pre-warning and problem indicator.

Mechanical Overload

If the pump cannot reach the set speed, its LED will **flash alternating yellow – green**. The pump will continue to work at reduced speed to e.g. run-in a new (hard!) tubing or to remove congestions from the piping system. If the drive still cannot reach the set speed after max. 5 min, the pump shuts down and the LED changes to **red**.

Thermal Overload

If the pump drive gets critically hot, the LED first **flashes yellow – red**.

Should the temperature rise further, the pump shuts down and the LED **flashes red**.

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Overload handling

Electrical Overload

If the mains voltage is too low and the pump cannot reach the set speed, its LED **flashes alternating yellow – green** (see ‚Mechanical Overload‘ above).

In case of mains voltage being less than nominal – 10 % but the drive has no problem to operate at the set speed, the LED **flashes green** to inform the user. The pump continues to run without restrictions.

⚠ **Before** the pump is restarted after any overload shut-down, the cause of the overload has to be determined (e.g. excessive differential pressure, congestion, cooling down of a high viscosity medium, wrong tubing wall thickness, instable mains voltage, environmental temperature, blocking of the rollerhead due to broken tubing).

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Status Indication

LED	Operation	Status class	Cause	Solution
green	Pump runs normally	normal		No fault
yellow	Pump stand-by	normal		No fault
flashes green-yellow	Adaptive speed ramp-up	Info	The set speed requires more torque than available. Thus, the actual speed is adaptively ramped up (duration: max. 5 min). <ol style="list-style-type: none"> 1 The backpressure is too high. 2 Under the application conditions, peak pressures occur during start-up. 3 A new (hard!) tube is being used. 4 The tubing had been fully squeezed for a long time and suffers deformation. 5 Tubing wall thickness is larger than 1.6 mm (1/16 inch). 6 A high viscosity medium is being pumped. 7 Mains voltage is less than nominal - 10%. 	Let the pump continue to run! If it shuts down with a red LED: check operation conditions.
flashes red-yellow	Temperature warning	Info	The internal temperature of the pump is becoming critical. <ol style="list-style-type: none"> 1 The fan is broken. 2 The environmental temperature is above 40°C (104°F). 	<ol style="list-style-type: none"> 1 Check if the fan expels air through the ventilation slots on the side. If it doesn't, contact your ISMATEC service center. 2 Provide for enough ventilation around the pump.
flashes green	Undervoltage	Info	Mains voltage is less than nominal - 10%.	Check the mains supply.
flashes yellow	Pumphead opened up while running	Fault	The pumphead was opened while the pump was running. <ol style="list-style-type: none"> 1 For safety reasons, the lid should be opened only when the pump is at rest. 2 The closed lid is not perfectly in position because the tubing is missing. 	<ol style="list-style-type: none"> 1 Close the lid and press the control knob once. The pump switches to 'stand-by'. 2 Insert a tube and close the lid. Press the control knob once to switch the pump to 'stand-by'.
flashes red	Overtemperature shut-down	Fault	Due to overheating, the drive shut down automatically. <ol style="list-style-type: none"> 1 The fan is broken. 2 The ventilation slots are blocked. 3 The environmental temperature is above 40°C. 	Press the control knob once. If the LED changes to yellow ('stand-by'), the pump has cooled down already. If not, wait for a couple of minutes. In any case, check if: <ol style="list-style-type: none"> 1 the fan expels air through the slots on the side. If it doesn't, contact you ISMATEC service center. 2 the ventilation slots on the side and in the bottom are blocked. Clean them if necessary. 3 the environmental temperature is above 40°C.
red	Overload shut-down	Fault	Maximum torque has been exceeded: <ol style="list-style-type: none"> 1 The drive is blocked. 2 An unsuitable tubing was inserted. 3 The tubing is not mounted correctly. 4 System pressure has changed and exceeds the specified maximum differential pressure. 5 Viscosity has gone up. 6 Mains voltage is less than nominal - 10%. 	<ol style="list-style-type: none"> 1 Make sure that the pumphead is not blocked. 2 Only suitable peristaltic tubing must be used (diameter, material, wall thickness - see page 22). 3 Ascertain that the tubing has been properly inserted (see page 13). 4 Verify that no excess backpressure or vacuum has build up at the outlet or inlet (e.g. by congestion, valves, filters, empty reservoirs without venting). 5 Inspect the medium for its temperature or other physical or chemical changes. 6 Check the mains supply. Press the control knob once to reset the pump. Its LED changes to stand-by
flashes yellow after power-on	Transmission of diagnostic data	normal	Right after power-up, the pump sends service information through the LED for service purposes.	No fault

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⚠ Flow rates per channel

In the tables listed subsequently the values indicated are only approximate and determined as follows:
ml/min, per channel, with water (22°C) and Tygon® tubing, without differential pressure

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Tubing materials wall thickness 1.6 mm

Tygon® LFL
Tygon® ST R-3603 / R-3607
PharMed® (Ismaprene)
Tygon® HC F-4040-A
Tygon® MH 2075
Tygon® SI Silicone 3350
Silicone Peroxid
C-Flex® 50-A
Norprene® A-60-G
Fluran HCA F-5500-A

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

Provided the tubing pump *Ecofluidic* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective *Ecofluidic* to your ISMATEC® representative. Please give information on defect, date of purchase, serial-no., and model.

Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask your ISMATEC® distributor for:

- parts lists
- replacement parts
- wiring diagrams
- spare brushes for motor

Please give information on defect, date of purchase, serial-no., and model.

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Cleaning

Materials

The dark plastic parts of the housing are made of PPS, the transparent colored ones of SAN (exceptions: power socket and housing feet). PPS has a high resistance against many chemicals. Nevertheless, we recommend not to use organic solvents for cleaning purposes but rather a common detergent.

Metal parts of the housing are made of stainless steel. Rollerhead, tubebed and knurled screw of the tubing fixation consist of anodized aluminum and stainless steel.

➔ The printed label with model and serial number has limited chemical resistance.

Protection rating of the housing is IP31. Therefore, it should be cleaned with a wet cloth only.

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
Cleaning

Disassembly of the pumphead

Switch off the pump and disconnect the power cable from mains.

- 1 The pumphead lid can be removed with an Allen key 3 mm. It suffices to loosen the Allen key screw with a few turns. After cleaning just put the lid back on its shaft, push down to the end and tighten the Allen key screw.
- 2 Similarly, the rollerhead can be removed with an Allen key 2.5 mm. Be cautious to maintain the greasing of the roller ball bearings when cleaning it. For assembly, just push the rollerhead onto its shaft down to the end and lock the Allen key screw.

The closure lock lever can simply be taken off its shaft and put back.

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Differential pressure

1) Possible differential pressure with appropriate tubing material, tubing with small i.d.'s may enable higher pressures

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< 70 db (A)

Materials

Housing made of PPS (GF40) and stainless steel, lid and control knob made of SAN

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Please note

We also recommend you to observe our **Warranty Terms as well as our Terms and Conditions of Sale**.

In case of any queries, please **contact your local ISMATEC® tmrepresentative**.

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Technical Specifications

Operating conditions

Temperature +41 to +104°F
Rel. humidity max. 80%
– not condensing, at normal laboratory conditions

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Ask for
the latest
ISMATEC
catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Laboratoriumstechnik GmbH