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#### Important information

##### ■ Safety

The installation of this panel-mounted pump must be carried out by a skilled person who is aware of the potential hazard involved in handling electrical components.

##### ■ Wiring

See figure on front page as well as the attached wiring diagram.

##### ■ Active cooling system

Depending on the installation and the operating time an active cooling system (ventilator) must be installed.

##### ■ Connecting cable

Shielded cables are to be used between control unit and drive. The shield must be connected to ground on both ends.

The maximum length of the cable checked and approved by Ismatec is 2 metres. Depending on the situation, longer cables can be used, which however may cause malfunction. For more information, please do not hesitate to contact us.

##### ■ CE conformity

Our standard pump models comply with the applicable EC directives.

However within the EC, our panel-mounted pump models must not be operated before the complete system or instrument complies with all the corresponding EC directives.

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Pins	Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0–5V	100 kΩ	ON*	OFF*	OFF*	OFF*	ON*
	0–10V	16.7 kΩ	OFF	ON	OFF	OFF	ON
	0–20mA	250 Ω	OFF	OFF	ON	OFF	ON
	4–20mA	250 Ω	OFF	OFF	OFF	ON	OFF

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## Analog interface

### Pin 1, (GND) Ground

Reference potential for all other inputs.

### Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

### Pin 4, direction

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

### Pin 4, direction

If this pin is left open, the pump runs in clockwise direction. Connecting pin 4 to the ground (pin1) makes the pump run in counter-clockwise direction.

### Pin 5, speed IN

External speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Input impedance and input range can be selected via a dip-switch inside the pump (see below).

### Pin 7, +36 V<sub>DC</sub>

About +36 V<sub>DC</sub> (unstabilised) are available the pump motor (max. load 1 A).

### Pin 9, speed OUT

The default setting is 0–11.2 V<sub>DC</sub>, proportionally to the motor speed 3.5–350 rpm. Alternatively a frequency range from 0 to 11.2 kHz is available. The input range can be selected via a dip-switch inside the pump (see below).

### Pin 10, +5V<sub>DC</sub>

This pin is connected to the stabilised +5V<sub>DC</sub> voltage supply of the pump.

### Pin 12, idle, overload

This pin is at +5V<sub>DC</sub> while the pump runs. When idle or in overload, the voltage changes to 0V<sub>DC</sub>.

➔ Switch off the mains voltage and let the device cool down for 2 minutes.

### Pin 13, internal speed

If pin 13 is connected with pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at pin 5.

## Settings of switch S1

\* Default setting

### Sliding switch S2

This switch affects Pin 9, speed OUT

Position A: 0–10 V<sub>DC</sub> (standard)  
Position B: 0–11.2 kHz

➔ Make sure that the pump is disconnected from the mains supply.

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