

Chemical Resistance Properties of Tubing

IMPORTANT: It is the user's responsibility to insure the suitability and safety of tubing for all intended uses. Laboratory, field or clinical tests must be conducted in accordance with applicable requirements in order to determine the safety and effectiveness for use of tubing in any particular application.

All ratings are based on room temperature (21°C; 73°F)

1 Excellent, 2 Good, 3 Fair, 4 Not Recommended Environment, % Conc. w-Water alc-Alcohol

	Tygon LFL	Tygon ST R-3603/R-3607	Pharmed Ismaprene	Tygon HC F-4040-A	Tygon MHLL	Tygon MHSL	Tygon SI Silikon 3350	Norprene A-60-G	Fluran HCA F-5500-A	CHEM-SURE	STA-PURE
							Platin				
Acetaldehyde	4	4	4	4	3	4	3	4	4	1	2
Acetamide, 67% in w	4	4	2	4	1	2	1	2	4	-	2
Acetate Solvents	4	4	2	4	4	4	4	2	4	-	-
Acetic Acid, 10% in w	1	1	1	1	1	2	1	1	4	1	2
Acetic Acid, 50-60% in w	1	1	2	1	1	1	1	2	4	1	2
Acetic Acid, Glacial, 100%	4	4	2	4	1	2	4	2	4	1	2
Acetic Anhydride	4	4	1	4	1	1	1	1	4	1	4
Acetone	4	4	4	4	2	3	3	4	4	1	4
Acetonitrile	4	4	2	4	2	2	4	2	4	-	-
Acetyl Bromide	4	4	3	4	4	3	4	3	4	-	-
Acetyl Chloride	4	4	3	4	4	3	4	3	4	1	-
Acetylene Gas	1	1	1	1	1	1	1	1	1	1	-
Acrylonitrile	4	4	2	4	2	2	4	2	4	1	4
Adipic Acid, 100% in alc	4	4	2	3	4	2	4	2	4	1	-
Air	1	1	1	1	1	1	1	1	1	1	1
Alcohols General	4	4	1	2	1	1	2	1	4	1	-
Aliphatic Hydrocarbons	4	4	4	2	4	4	4	4	2	-	-
Allyl Alcohol	4	4	3	1	1	1	4	3	1	-	-
Alum, 5% in w	1	1	1	1	1	1	1	1	1	-	-
Aluminum Chloride, 53% in w	1	1	1	1	1	1	1	1	1	-	-
Aluminum Hydroxide, 2% in w	1	1	1	1	1	1	1	1	1	-	-
Aluminum Sulfate, 50% in w	1	1	1	1	1	1	1	1	1	-	2
Aluminum Salts	1	1	1	1	1	1	1	1	1	-	2
Amines	4	4	3	4	4	4	4	3	4	2	-
Ammonia Gas	1	1	1	1	1	1	4	1	4	1	2
Ammonia, Anhydrous Liquid	2	2	2	2	2	2	4	2	4	-	-
Ammonium Acetate, 45% in w	1	1	1	1	1	1	1	1	4	1	-
Ammonium Carbonate, 20% in w	1	1	1	1	1	1	1	1	1	1	-
Ammonium Hydroxide, 5-10% in w	1	1	1	2	1	1	4	1	4	1	2
Ammonium Hydroxide, 30% in w	1	1	1	3	1	1	4	1	4	1	2
Ammonium Persulfate, 30% in w	1	1	1	1	1	1	1	1	1	-	-
Ammonium Salts	1	1	1	1	1	1	1	1	1	-	-
Ammonium Sulfate, 30% in w	1	1	1	1	1	1	1	1	1	1	2
Amyl Acetate	4	4	2	4	4	4	4	2	4	2	4
Amyl Alcohol	4	4	4	1	1	1	4	4	1	1	4
Amyl Chloride	4	4	3	4	4	4	4	3	4	-	4
Aniline	4	4	3	4	4	4	4	3	4	1	2
Aniline Hydrochloride	4	4	3	4	4	4	4	3	4	1	4
Antimony Salts	1	1	1	1	1	1	1	1	1	-	-
Aqua Regia	4	4	4	4	1	1	4	4	4	-	-
Aromatic Hydrocarbons	4	4	4	4	4	4	4	4	4	-	4
Arsenic Acid, 20% in w	1	1	3	1	1	1	3	3	1	1	-
Arsenic Salts	1	1	1	1	1	1	1	1	1	-	-
ASTM Reference No. 1 Oil	4	4	3	1	4	4	1	3	1	-	-
ASTM Reference No. 2 Oil	4	4	4	1	4	4	2	4	1	-	-
ASTM Reference No. 3 Oil	4	4	4	1	4	4	4	4	1	-	-
Barium Carbonate, 1% in w	1	1	1	1	1	1	1	1	1	-	-
Barium Hydroxide, 5% in w	1	1	1	1	1	1	1	1	1	1	2
Beer	1	1	1	1	1	1	1	1	1	1	2
Benzaldehyde	4	4	4	4	3	3	3	4	4	1	4
Benzene	4	4	4	4	4	4	4	4	4	1	4
Benzenesulfonic Acid	4	4	4	4	4	4	4	4	4	1	-
Benzoic Acid	4	4	2	3	4	4	4	2	4	1	2
Benzyl Alcohol	4	4	1	4	1	1	1	1	1	-	-
Bleach Liquor, 22% in w	1	1	1	1	1	1	4	1	1	1	-
Borax, 6% in w	1	1	1	1	1	1	1	1	1	-	2
Boric Acid, 4% in w	1	1	1	1	1	1	1	1	1	1	2
Bromine, Anhydrous Liquid	4	4	4	4	4	4	4	4	4	-	4
Butadiene	1	1	1	1	2	2	1	1	1	2	-
Butane	1	1	1	1	2	2	1	1	1	2	4
Butyl Acetate	4	4	2	4	4	4	4	2	4	2	4
Butyl Alcohol	4	4	4	1	1	1	4	4	1	1	2
Butyric Acid	4	4	2	3	4	4	4	2	4	1	-
Calcium Carbonate, 25% in dilute acids	1	1	1	1	1	1	1	1	1	1	-
Calcium Chloride, 30% in w	1	1	1	1	1	1	1	1	1	-	2
Calcium Hydroxide, 10% in glycerol	1	1	1	4	1	1	1	1	1	1	2
Calcium Hypochlorite, 20% in w	1	1	1	1	1	1	4	1	1	1	2
Calcium Nitrate, 55% in w	1	1	1	1	1	1	1	1	1	1	2
Calcium Salts	1	1	1	1	1	1	1	1	1	-	2
Calcium Sulfate, 0.2% in w	1	1	1	1	1	1	1	1	1	1	-
Carbon Dioxide, Wet/Dry	1	1	1	1	1	1	1	1	1	1	2
Carbon Disulfide	4	4	4	4	4	4	4	4	2	-	-
Carbon Monoxide	1	1	1	1	1	1	1	1	1	1	2
Carbon Tetrachloride	4	4	4	4	4	4	4	4	4	2	4
Carbonic Acid	1	1	1	1	1	1	1	1	1	1	2
Castor Oil	4	4	3	1	2	2	1	3	1	1	2
Cellosolve	4	4	3	3	4	4	4	3	4	-	-
Cellosolve Acetate	4	4	3	3	4	4	4	3	4	-	-
Chlorine, Dry Gas	1	1	3	1	3	3	4	3	1	-	4

	Tygon LFL	Tygon ST R-3603/R-3607	Pharmed Ismaprene	Tygon HC F-4040-A	Tygon MHLL	Tygon MHSL	Tygon S1 Silikon 3350 Platin	Norprene A-60-G	Fluran HCA F-5500-A	CHEM-SURE	STA-PURE
Chlorine, Wet Gas	2	2	4	1	3	3	4	4	1	-	4
Chloroacetic Acid, 20% in w	1	1	2	4	1	1	1	2	4	2	-
Chlorobenzene, Mono, Di, Tri	4	4	4	4	4	4	4	4	4	1	4
Chloroform	4	4	3	4	4	4	4	3	4	2	4
Chlorosulfonic Acid	4	4	4	4	4	4	4	4	4	1	4
Chromic Acid, 10-20% in w	2	2	1	3	1	1	4	1	1	-	4
Chromic Acid, 50% in w	3	3	3	4	2	2	4	3	1	-	4
Citric Acid, 10-20% in w	1	1	1	1	1	1	1	1	1	1	2
Coconut Oil	4	4	3	1	2	2	1	3	1	1	2
Corn Syrup	1	1	1	1	1	1	1	1	1	-	-
Cottonseed Oil	4	4	3	1	2	3	1	3	1	1	-
Cresol (m, o, or p)	4	4	4	3	1	1	2	4	1	1	4
Cresylic Acid	4	4	2	3	4	4	4	2	4	-	-
Cupric Chloride, 40% in w	1	1	1	1	1	1	1	1	1	-	-
Cupric Nitrate, 70% in w	1	1	1	1	1	1	1	1	1	-	-
Cupric Sulfate, 13% in w	1	1	1	1	1	1	1	1	1	-	-
Cyclohexane	4	4	4	3	4	4	4	4	1	2	4
Cyclohexanone	4	4	4	4	3	3	4	4	4	-	-
Detergent Solutions	1	1	2	1	1	1	1	2	1	1	2
Dibutyl Phthalate	3	3	1	3	1	1	1	1	1	1	2
Diesel Fuel	4	4	4	2	4	4	4	4	1	1	-
Diethylamine	1	1	1	3	1	1	4	1	4	2	2
Diethylene Glycol	1	1	1	1	1	1	1	1	1	1	2
Dimethylformamide	4	4	2	4	1	1	1	2	4	1	2
Dimethylsulfoxide	4	4	1	4	2	2	3	1	4	-	-
Diethyl Phthalate	3	3	1	3	1	1	1	1	1	1	4
Dioxane	4	4	4	4	4	4	4	4	4	-	2
Ether	4	4	3	3	4	4	4	3	4	-	2
Ethyl Acetate	4	4	2	4	4	4	4	2	4	1	4
Ethyl Alcohol	4	4	3	2	1	1	3	3	4	1	4
Ethyl Benzoate	4	4	4	4	4	4	4	4	4	-	-
Ethyl Chloride	4	4	3	4	4	4	4	3	4	-	4
Ethyl Ether	4	4	3	3	4	4	4	3	4	2	-
Ethylene Bromide	4	4	4	4	3	3	1	4	1	-	-
Ethylene Chlorohydrin	4	4	1	2	1	1	2	1	4	-	4
Ethylene Dichloride	4	4	3	4	4	4	4	3	4	-	-
Ethylene Glycol	1	1	1	1	1	1	1	1	1	1	2
Ethylene Oxide	1	1	1	1	1	1	1	1	1	2	4
Fatty Acids	4	4	3	2	3	3	2	3	1	1	-
Ferric Chloride, 43% in w	1	1	1	1	1	1	1	1	1	-	2
Ferric Nitrate, 60% in w	1	1	1	1	1	1	1	1	1	-	-
Ferric Sulfate, 5% in w	1	1	1	1	1	1	1	1	1	1	2
Ferrous Chloride, 40% in w	1	1	1	1	1	1	1	1	1	-	-
Ferrous Sulfate, 5% in w	1	1	1	1	1	1	1	1	1	1	-
Fluoboric Acid, 48% in w	1	1	4	4	1	1	4	4	1	-	-
Fluorine Gas	4	4	4	4	4	4	4	4	2	4	-
Fluosilicic Acid, 25% in w	1	1	3	1	1	1	3	3	1	-	-
Formaldehyde, 37% in w	4	4	4	4	3	3	3	4	4	1	2
Formic Acid, 25% in w	1	1	1	3	1	1	1	1	1	1	2
Formic Acid, 40-50% in w	2	2	2	4	1	1	1	2	3	1	2
Formic Acid, 98% in w	2	2	2	4	1	1	1	2	4	1	2
Freon 11	1	1	1	1	1	1	1	1	2	4	-
Freon 12	1	1	1	1	1	1	1	1	2	4	-
Freon 22	1	1	1	1	1	1	1	1	2	4	-
Fruit Juice	1	1	1	1	1	1	1	1	1	1	-
Fuel Oil	4	4	4	2	4	4	4	4	1	1	4
Furfural	4	4	4	4	4	4	4	4	4	2	-
Gallic Acid, 17% in acetone	4	4	2	3	4	4	4	2	4	-	-
Gasoline, Automotive	4	4	4	2	4	4	4	4	1	2	4
Gelatin	1	1	1	1	1	1	1	1	1	1	2
Glucose, 50% in w	1	1	1	1	1	1	1	1	1	1	2
Glycerin	1	1	1	1	1	1	1	1	1	1	2
Glycolic Acid, 70% in w	1	1	2	1	1	1	1	2	4	-	-
Heptane	4	4	4	2	4	4	4	4	2	2	-
Hexane	4	4	4	2	4	4	4	4	2	-	4
Hydrazine	4	4	3	4	4	4	4	3	4	1	4
Hydrobromic Acid, 20-50% in w	1	1	4	1	1	1	4	4	1	1	4
Hydrobromic Acid, 100% in w	1	1	4	4	1	1	4	4	1	1	4
Hydrochloric Acid, 10% in w	1	1	1	1	1	1	1	1	1	1	2
Hydrochloric Acid, 37% in w	1	1	2	4	1	1	4	2	2	1	4
Hydrocyanic Acid	1	1	1	1	1	1	1	1	4	1	2
Hydrofluoric Acid, 10% in w	1	1	4	1	1	1	4	4	1	-	2
Hydrofluoric Acid, 25% in w	1	1	4	4	1	1	4	4	1	-	2
Hydrofluoric Acid, 40-48% in w	3	2	4	4	1	1	4	4	1	-	2
Hydrogen Gas	1	1	1	1	1	1	1	1	1	1	4
Hydrogen Peroxide, 3% in w	1	1	1	1	1	1	1	1	1	1	2
Hydrogen Peroxide, 10% in w	1	1	1	1	1	1	1	1	1	1	2
Hydrogen Peroxide, 30% in w	1	1	1	4	1	1	1	1	1	1	2
Hydrogen Peroxide, 90% in w	4	4	2	4	2	2	3	2	1	1	2
Hydrogen Sulfide	1	1	1	1	1	1	1	1	1	1	4
Hydroquinone, 7% in w	1	1	2	1	1	1	2	2	1	-	-
Hypochlorous Acid, 25% in w	1	1	1	1	1	1	1	1	1	1	-
Iodine, 50 ppm in w	1	1	1	1	1	1	1	1	1	-	-
Isobutyl Alcohol	4	4	3	1	1	1	4	3	1	-	2
Isooctane	4	4	4	2	4	4	4	4	1	-	4
Isopropyl Acetate	4	4	2	4	4	4	4	2	4	2	-
Isopropyl Alcohol	4	4	3	1	1	1	4	3	1	1	2

	Tygon LFL	Tygon ST R-3603/R-3607	Pharmed Ismaprene	Tygon HC F-4040-A	Tygon MHLL	Tygon MHSL	Tygon S1 Silikon 3350 Platin	Norprene A-60-G	Fluran HCA F-5500-A	CHEM-SURE	STA-PURE
Isopropyl Ether	4	4	3	3	4	4	4	3	4	2	-
Jet Fuel, JP8	4	4	4	2	4	4	4	4	1	1	-
Kerosene	4	4	4	2	4	4	4	4	1	1	4
Ketones	4	4	4	4	3	3	4	4	4	-	4
Lacquer Solvents	4	4	2	4	4	4	4	2	4	1	-
Lactic Acid, 3-10% in w	1	1	1	1	1	1	1	1	4	1	2
Lactic Acid, 85% in w	4	4	2	4	1	1	4	2	4	1	2
Lard, Animal Fat	4	4	3	1	2	2	1	3	1	1	2
Lead Acetate, 35% in w	1	1	1	1	1	1	1	1	1	-	4
Lead Salts	1	1	1	1	1	1	1	1	1	-	-
Lemon Oil	4	4	4	3	4	4	4	4	1	-	-
Limonene-D	4	4	4	3	4	4	4	4	1	-	-
Linoleic Acid	4	4	3	2	3	3	2	3	1	1	-
Linseed Oil	4	4	3	1	2	2	1	3	1	1	2
Lubricating Oils, Petroleum	4	4	4	1	4	4	2	4	1	1	4
Magnesium Carbonate, 1% in w	1	1	1	1	1	1	1	1	1	-	-
Magnesium Chloride, 35% in w	1	1	1	1	1	1	1	1	1	1	2
Magnesium Hydroxide, 10% in dil. acid	1	1	1	1	1	1	1	1	1	1	-
Magnesium Nitrate, 50% in w	1	1	1	1	1	1	1	1	1	-	-
Magnesium Sulfate, 25% in w	1	1	1	1	1	1	1	1	1	1	2
Maleic Acid, 30% in w	4	4	3	2	3	3	2	3	1	1	-
Malic Acid, 36% in w	1	1	1	1	1	1	1	1	4	1	-
Manganese Salts	1	1	1	1	1	1	1	1	1	-	2
Mercuric Chloride, 6% in w	1	1	1	1	1	1	1	1	1	1	-
Mercuric Cyanide, 8% in w	1	1	1	1	1	1	1	1	1	-	-
Mercury	1	1	1	1	1	1	1	1	1	1	2
Mercury Salts	1	1	1	1	1	1	1	1	1	-	-
Methane Gas	1	1	1	1	1	1	1	1	1	2	4
Methanol	4	4	4	2	1	4	2	1	4	1	2
Methyl Acetate	4	4	2	4	4	4	4	2	4	2	-
Methyl Bromide	4	4	3	4	4	4	4	3	4	2	2
Methyl Chloride	4	4	3	4	4	4	4	3	4	2	4
Methyl Ethyl Ketone	4	4	4	4	3	3	4	4	4	2	4
Methyl Isobutyl Ketone	4	4	4	4	3	3	4	4	4	2	4
Methylene Chloride	4	4	3	4	4	4	4	3	4	2	4
Methyl Methacrylate	4	4	4	4	4	4	4	4	3	1	2
Milk	1	1	1	1	1	1	1	1	1	1	2
Mineral Oil	3	3	4	1	4	4	4	4	1	1	2
Mineral Spirits	4	4	4	2	4	4	4	4	1	-	-
Molasses	1	1	1	1	1	1	1	1	1	1	-
Monoethanolamine	4	4	3	4	4	4	4	3	4	-	-
Motor Oil	4	4	4	1	4	4	4	4	1	-	-
Naphtha	4	4	4	2	4	4	4	4	1	2	4
Naphthalene	4	4	4	2	4	4	4	4	1	1	4
Natural Gas	1	1	1	1	1	1	1	1	1	2	2
Nickel Chloride, 40% in w	1	1	1	1	1	1	1	1	1	1	2
Nickel Nitrate, 75% in w	1	1	1	1	1	1	1	1	1	-	-
Nickel Salts	1	1	1	1	1	1	1	1	1	-	2
Nickel Sulfate, 25% in w	1	1	1	1	1	1	1	1	1	1	-
Nitric Acid, 10% in w	1	1	1	4	1	1	3	1	1	1	4
Nitric Acid, 35% in w	1	1	1	4	1	1	4	1	3	1	4
Nitric Acid, 68-71% in w	4	4	4	4	1	1	4	4	4	-	4
Nitrobenzene	4	4	4	4	4	4	4	4	4	1	4
Nitromethane	4	4	4	4	4	4	4	4	4	1	-
Nitrous Acid, 10% in w	1	1	1	3	1	1	2	1	1	-	-
Nitrous Oxide	1	1	1	1	1	1	1	1	1	-	-
Oils, Animal	4	4	3	1	2	2	1	3	1	-	2
Oils, essential	4	4	4	3	4	4	4	4	1	-	-
Oils, Hydraulic (Phosphate Ester)	4	4	4	1	4	4	4	4	1	-	-
Oils, Hydrocarbon	4	4	4	1	4	4	2	4	1	-	-
Oils, Vegetable	4	4	3	1	2	2	1	3	1	1	2
Oleic Acid	4	4	3	2	3	4	2	3	1	1	4
Oleum, 25% in w	1	1	1	2	1	1	2	1	1	-	-
Ortho Dichlorobenzene	4	4	4	4	4	4	4	4	4	-	-
Oxalic Acid, 12% in w	2	2	2	4	1	1	1	2	4	1	2
Oxygen	1	1	1	1	1	1	1	1	1	1	2
Ozone, 300pphm	1	1	1	1	1	1	1	1	1	1	2
Palmitic Acid, 100% in ether	4	4	3	2	3	3	2	3	1	1	-
Paraffins	4	4	4	2	4	4	4	4	2	1	4
Perchloric Acid, 67% in w	1	3	1	4	1	1	4	1	1	1	4
Perchloroethylene	3	4	3	4	4	4	4	3	4	2	2
Phenol, 5-10% in w	4	2	1	1	1	1	4	1	1	1	4
Phenol, 91% in w	4	4	1	3	1	1	4	1	1	-	4
Phosphoric Acid, <10% in w	1	1	1	1	1	1	3	1	1	1	4
Phosphoric Acid, 25% in w	1	1	1	1	1	1	4	1	1	1	4
Phosphoric Acid, 85% in w	2	1	1	4	1	1	4	1	1	1	4
Phosphorous Trichloride Acid	2	1	2	4	1	1	4	2	2	-	-
Photographic Solutions	2	1	2	1	1	1	2	2	1	-	-
Phthalic Acid, 9% in alc	4	4	1	3	1	1	2	1	1	-	2
Phthalic Anhydride, 9% in alc	2	4	1	4	1	1	1	1	4	1	-
Picric Acid, 1% in w	4	1	4	1	1	1	4	4	1	1	4
Plating Solutions	2	1	1	4	1	1	4	1	1	-	-
Potassium Carbonate, 55% in w	1	1	1	1	1	1	1	1	1	1	2
Potassium Cyanide, 33% in w	1	1	1	1	1	1	1	1	1	1	2
Potassium Dichromate, 5% in w	1	1	1	1	1	1	1	1	1	1	2
Potassium Hydroxide, <10% in w	1	1	1	4	1	1	1	1	1	1	4
Potassium Hypochlorite, 70% in w	1	1	1	1	1	1	1	1	1	-	-

	Tygon LFL	Tygon ST R-3603/R-3607	Pharmed Ismaprene	Tygon HC F-4040-A	Tygon MHLL	Tygon MHSL	Tygon S1 Silikon 3350 Platin	Norprene A-60-G	Fluran HCA F-5500-A	CHEM-SURE	STA-PURE
Potassium Iodide, 56% in w	1	1	1	1	1	1	1	1	1	-	2
Potassium Permanganate, 6% in w	1	1	1	1	1	1	1	1	1	1	-
Potassium Salts	1	1	1	1	1	1	1	1	1	-	2
Propane Gas	1	1	1	1	1	1	1	1	1	2	4
Propyl Alcohol (Propanol)	4	4	3	1	1	1	4	3	1	-	-
Propylene Glycol	1	1	1	1	1	1	1	1	1	-	-
Propylene Oxide	1	1	1	1	1	1	1	1	1	2	-
Pyridine	1	4	3	4	3	3	4	3	4	1	4
Salicylic Acid, 1% in w	1	1	1	1	1	1	1	1	4	-	-
Silicone Oils	1	2	3	1	1	2	4	3	1	-	4
Silver Nitrate, 55% in w	1	1	1	1	1	1	1	1	1	1	2
Skydrol 500A	4	4	4	1	4	4	4	4	1	-	-
Soap Solutions	1	1	2	1	1	1	1	2	1	1	2
Sodium Acetate, 55% in w	1	1	1	1	1	1	1	1	1	1	-
Sodium Benzoate, 22% in w	1	1	1	1	1	1	1	1	1	-	-
Sodium Bicarbonate, 7% in w	1	1	1	1	1	1	1	1	1	1	2
Sodium Carbonate, 7% in w	1	1	1	1	1	1	1	1	1	1	-
Sodium Chlorate, 45% in w	1	1	1	1	1	1	1	1	1	-	-
Sodium Chloride, 20% in w	1	1	1	1	1	1	1	1	1	1	2
Sodium Cyanide, 30% in w	1	1	1	4	1	1	1	1	1	1	2
Sodium Fluoride, 3% in w	1	1	1	1	1	1	1	1	1	-	-
Sodium Hydroxide, 10-15% in w	1	1	1	4	1	1	1	1	1	1	2
Sodium Hydroxide, 30-40% in w	1	3	1	4	1	1	1	1	1	-	2
Sodium Hypochlorite, 5.5% in w	1	1	1	1	1	1	4	1	1	1	2
Sodium Hypochlorite, 12.2% in w	1	1	1	1	1	1	4	1	1	1	2
Sodium Nitrate, 3.5% in w	1	1	1	1	1	1	1	1	1	1	4
Sodium Salts	1	1	1	1	1	1	1	1	1	-	-
Sodium Sulfates, 3.6% in w	1	1	1	1	1	1	1	1	1	1	2
Sodium Sulfide, 13% in w	1	1	1	1	1	1	1	1	1	1	2
Sodium Sulfite, 10% in w	1	1	1	1	1	1	1	1	1	-	-
Stannic Chloride, 50% in w	1	1	1	1	1	1	1	1	1	-	2
Stannous Chloride, 45% in w	1	1	1	1	1	1	1	1	1	-	2
Stearic Acid, 5% in alc	1	4	3	2	3	3	2	3	1	1	2
Styrene Monomer	4	4	4	4	4	4	4	4	3	1	-
Sulfur Chloride	1	4	4	4	1	1	4	4	1	1	4
Sulfur Dioxide, Dry Gas	1	1	1	2	1	1	1	1	1	1	2
Sulfur Dioxide, Wet gas	1	1	1	2	1	1	1	1	1	1	2
Sulfur Trioxide, Wet	1	2	2	4	2	2	2	2	2	-	2
Sulfuric Acid, 10% in w	1	1	1	1	1	1	1	1	1	1	4
Sulfuric Acid, 30% in w	1	1	1	2	1	1	2	1	1	1	4
Sulfuric Acid, 95-98% in w	4	4	4	4	1	4	4	4	1	1	4
Sulfurous Acid	1	1	1	1	1	1	1	1	1	1	4
Tannic Acid, 75% in w	1	2	2	4	1	1	1	2	4	-	2
Tartaric Acid, 56% in w	1	1	1	1	1	1	1	1	1	1	2
Tetrahydrofuran	4	4	4	4	4	4	4	4	4	2	-
Thionyl Chloride	1	1	1	2	1	1	1	1	1	-	4
Tin Salts	1	1	1	1	1	1	1	1	1	-	-
Titanium Salts	1	1	1	1	1	1	1	1	1	-	-
Toluene	4	4	4	4	4	4	4	4	3	1	4
Trichloroacetic Acid, 90% in w	1	1	2	4	1	1	1	2	4	-	-
Trichloroethane	4	4	3	4	4	4	4	3	4	2	4
Triethanolamine	4	4	3	4	4	4	4	3	4	1	-
Trichloroethylene	4	4	4	4	4	4	4	4	4	2	4
Trichloropropane	4	4	3	4	4	4	4	3	4	-	-
Tricresyl Phosphate	2	3	1	3	1	1	1	1	1	1	4
Trisodium Phosphate	1	1	1	1	1	1	1	1	1	-	-
Turpentine	4	4	4	2	4	4	4	4	1	1	4
Urea, 20% in w	1	1	1	1	1	1	1	1	1	1	2
Uric Acid	1	1	1	3	1	1	1	1	3	-	-
Vinegar	1	1	1	1	1	1	1	1	4	1	2
Vinyl Acetate	4	4	2	4	4	4	4	2	4	2	-
Water, Deionized	1	1	1	1	1	1	1	1	1	1	1
Water, Distilled	1	1	1	1	1	1	1	1	1	1	1
Xylene	4	4	4	4	4	4	4	4	3	1	4
Zinc Chloride, 80% in w	1	1	1	1	1	1	1	1	1	1	-
Zinc Salts	1	1	1	1	1	1	1	1	1	-	-