


Control unit BVP-E Z / -Remote


for DC motor with encoder

OEM




Control unit BVP-E Z Remote
(Remote = without control panel)


Operating panel BVP-E Z



Front view (mm)



Mounting cut-out (mm)



IMPORTANT

Active cooling system
Adequate ventilation (active cooling) must be provided depending on the installation conditions and the operating time.
The buyer and user himself is responsible for sufficient ventilation of ISMATEC® panel-mounted pumps.

Technical specifications

Mains connection	230 V / 50/60 Hz or 115 V / 50/60 Hz (please indicate on your order)
Power consumption	max. 150 W
Weight	2.3 kg

Interfaces

■ Analog input	
Speed control	0-5 V or 0-10 V, 0-20 mA or 4-20 mA
■ Analog output	
Speed output	0-10 V _{dc} or 0-12 kHz
■ Digital input (TTL-level)	
Start/Stop	

Operating panel

- Start/Stop switch
- Switch for MAX/SLOW
- 3-digit potentiometer for digital speed setting in 0.1%-steps

Ordering information

Model	BVP-E Z	BVP-E Z Remote
	IS 2367	on request
Included with delivery: data sheet for installation and electrical connections		

Data sheet
08.11
30.09.05