


# Control unit BVP-E Z / -Remote


for DC motor with encoder

OEM




Control unit BVP-E Z Remote  
(Remote = without control panel)


Operating panel BVP-E Z



Front view (mm)



Mounting cut-out (mm)



**⚠ IMPORTANT**

**Active cooling system**  
Adequate ventilation (active cooling) must be provided depending on the installation conditions and the operating time.  
The buyer and user himself is responsible for sufficient ventilation of ISMATEC® panel-mounted pumps.

### Technical specifications

Mains connection	230 V / 50/60 Hz or 115 V / 50/60 Hz (please indicate on your order)
Power consumption	max. 150 W
Weight	2.3 kg

### Interfaces

■ Analog input	
Speed control	0-5 V or 0-10 V, 0-20 mA or 4-20 mA
■ Analog output	
Speed output	0-10 V <sub>dc</sub> or 0-12 kHz
■ Digital input (TTL-level)	
Start/Stop	

### Operating panel

- Start/Stop switch
- Switch for MAX/SLOW
- 3-digit potentiometer for digital speed setting in 0.1%-steps

### Ordering information

Model	BVP-E Z	BVP-E Z Remote
	IS 2367	on request
Included with delivery: data sheet for installation and electrical connections		

Data sheet  
08.11  
30.09.05